

PRACTICAL
MEDICINAL CHEMISTRY

-I-



JERASH UNIVERSITY
FACULTY OF PHARMACY
PRACTICAL MEDICINAL CHEMISTRY I
COURSE OUTLINE

LAB I:- Assay of Ammonium Chloride.

LAB II:- Assay of Ferrous Sulphate.

LAB III:- Assay of Promethazine.HCl Syrup.

LAB IV:- Assay of Mefenamic acid Tablets.

LAB V:- Assay of Ibuprofen Tablets.

LAB VI:- Assay of Aspirin Tablets and Salicylic acid in Aspirin.

LAB VII:- Assay of Frusemide Tablets.

LAB VIII:- Assay of Nalidixic acid Tablets.

LAB IX:- Assay of Dimenhydrinate Tablets.

LAB X:- Assay of Chlorpropamide Tablets.

Reports:- 10 marks

Quizes:- 10 marks

Mid-Term exam:- 30 marks

Final-Term exam:- 50 marks

JERASH UNIVERSITY
FACULTY OF PHARMACY
MEDICINAL CHEMISTRY I
LAB.I

ASSAY AND IDENTIFICATION OF AMMONIUM CHLORIDES

PART I-IDENTIFICATION TEST FOR CHLORIDE

ION:-1- Make solution of 0.1 gm Ammonium chloride in water 2- Acidify with 2M Nitric acid 3- Add 0.4 ml of silver nitrate solution 4- Shake and allow to stand, see the result and explain.

PART II-ACIDITY OR ALKALINITY OF AMMONIUM CHLORIDE:-

1-Dissolve 1 gm of Ammonium chloride in 10 ml of distilled water
2-Add 0.05 ml of Methyl red
3- Add 0.5 ml of 0.01M HCl, note the result
4- Add 1 ml of 0.01M NaOH, note the result
5- Write your conclusion about acidity or alkalinity of your solution.

PART III- ASSAY OF AMMONIUM CHLORIDE:-

1- Dissolve 1 gm of Ammonium Chloride in 20 ml of water
2- Add a mixture of 5 ml of Formaldehyde previously neutralized to Phenolphthalein solution and 20 ml of water
3- Allow to stand for 1-2 mins
4- Titrate slowly with 1M NaOH, using 0.2 ml of Phenolphthalein solution as an indicator.
The official standard of Ammonium chloride raw material is 99.0-100.5% according to the BP.

JERASH UNIVERSITY
FACULTY OF PHARMACY
MEDICINAL CHEMISTRY I
LAB.II

ASSAY AND IDENTIFICATION OF $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$

PART I- IDENTIFICATION OF IRON:-

- 1-Take about 10 mg of the salt. Dissolve in 1 ml of water
- 2- Add 1 ml of Potassium hexacyano ferrate(III) solution
- 3- Note the result. Add 2M HCl and see if any changes occur.

PART II- IDENTIFICATION OF SULPHATE:-1-

- Take about 45 mg of the salt. Dissolve in 5 ml of water
- 2- Add 1 ml of 2M HCl and 1 ml of Barium Chloride solution.
- Note the result.

PART III- ASSAY OF $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$:-

- 1- Dissolve 2.5 gm of sodium hydrogen carbonate in a mixture of 150 ml water and 10 ml Sulphuric acid
- 2- When effervescence ceases, add 0.5 gm of ferrous sulphate
- 3- Shake gently to dissolve and titrate with 0.1M Ammonium cerium (IV) nitrate using 0.1 ml of Ferroin solution as indicator.

JERASH UNIVERSITY
FACULTY OF PHARMACY
MEDICINAL CHEMISTRY I
LAB.III

ASSAY OF PROMETHAZINE.HCL SYRUP

Procedure:-

- 1- Take 10 ml of the syrup to be assayed volumetrically.
- 2- Add 5 ml of water and shake for 5 min.
- 3- Add 1.5 ml of 1.2N NaOH and shake for 5 min.
- 4- Extract with 2 successive quantities of ether, 7.5 ml each
- 5- Filter the organic layer then evaporate to dryness
- 6- Add 2.5 ml of alcohol to the residue and shake well to dissolve
- 7- Add 5 ml of 0.1 HCl volumetrically
- 8- Titrate the excess acid with 0.1N NaOH using methyl red indicator
- 9- Carry out the blank titration.
- 10- Calculate the amount of Promethazine.HCl in your sample.

Content of Promethazine.HCl should be 90-110% of the prescribed amount.

JERASH UNIVERSITY
FACULTY OF PHARMACY
MEDICINAL CHEMISTRY I
LAB.IV

ASSAY AND IDENTIFICATION OF MEFENAMIC ACID CAPSULES

PART I- IDENTIFICATION:-

I- Dissolve 25 mg of Mefenamic acid in 15 ml of Chloroform and examine under UV light (365 nm), the solution exhibits a strong greenish-yellow fluorescence, carefully add 0.5 ml of a saturated solution of trichloroacetic acid drop wise and examine under UV (365 nm) light, no fluorescence is seen

II- Dissolve 5 mg in 20 ml sulfuric acid and add 0.05 ml of 0.0167M potassium dichromate, an intense blue color is produced immediately, and it fades rapidly to brownish-green.

ASSAY OF MEFENAMIC ACID:-

1- Dissolve a quantity containing 250 mg of the drug in 50 ml of warm absolute ethanol previously neutralized to phenol red solution

2- Titrate with 0.1M sodium hydroxide using phenol red solution as indicator.

Each 1 ml of 0.1M sodium hydroxide is equivalent to 24.13 mg of $C_{15}H_{15}NO_2$

JERASH UNIVERSITY
FACULTY OF PHARMACY I
LAB.V

ASSAY OF IBUPROFEN TABLETS

Procrdure:-

- 1- Weigh and powder 20 tablets
- 2- Extract a quantity of tablets equivalent to 0.4 gm Ibuprofen with 20 ml of Chloroform for 15 min.
- 3- Filter under reduced pressure
- 4- Wash the residue with two 5 ml quantities of Chloroform
- 5- Evaporate the combined filtrate just to dryness in a current of air
- 6- Dissolve the residue in 80 ml of ethanol (previously neutralized to phenolphthalein indicator)
- 7- Titrate with 0.1M NaOH using 0.2 ml of phenolphthalein as an indicator
- 8- Repeat the same procedure with out the substance
- 9- The difference between the two titrations represents the amount of NaOH reacted with the drug.

Each 1 ml of 0.1M NaOH is equivalent to 20.63 mg Ibuprofen.

FACULTY OF PHARMACY
MEDICINAL CHEMISTRY I
LAB.VI

ASSAY OF ASPIRIN TABLETS

Procedure:-

- 1- Weigh and powder 20 tablets
- 2- To a quantity of the powder containing 0.3 gm of Aspirin add 30 ml of 0.5M sodium hydroxide VS
- 3- Boil gently for 10 min.
- 4- Titrate the excess of alkali with 0.5M HCl VS, using phenol red solution as indicator
- 5- Repeat the operation with out the substance being examined
- 6- The difference between the titration represents the amount of sodium hydroxide required.

Each ml of 0.5M sodium hydroxide VS is equivalent to 45.04 mg of $C_9H_8O_4$.

ASSAY OF SALICYLIC ACID IN ASPIRIN TABLETS

- 1- Shake a quantity of powdered tablets containing 0.2 gm of Aspirin in 4 ml of ethanol and dilute to 100 ml at a temperature not exceeding $10^{\circ}C$
- 2- Filter immediately and transfer 50 ml of the filtrate to Nessler cylinder
- 3- Add 1 ml of freshly prepared Iron III sulphate R1. Mix and allow to stand for 1 min
- 4- Any violet color produced is not more than that obtained by adding 1 ml of freshly prepared Ammonium iron III sulphate solution R1 to a mixture of 3 ml of freshly prepared 0.01% W/V solution of salicylic acid, 2 ml ethanol 96 % and sufficient water to produce 50 ml in a second Nessler tube.

JERASH UNIVERSITY
FACULTY OF PHARMACY
MEDICINAL CHEMISTRY I
LAB.VII

ASSAY OF FRUSEMIDE TABLETS

Procedure:-

- 1- Shake a quantity of the powder containing 0.04 of Frusemide with 60 ml of 0.1M NaOH for 10 min
- 2- Add sufficient 0.1M NaOH to produce 100 ml and filter.
- 3- Dilute 1 ml to 50 ml with 0.1M NaOH
- 4- Measure the absorbance at the maximum at 271 nm

Calculate the content of $C_{12}H_{11}ClN_2O_5S$ taking 580 as the value of A (1% ,1cm) at 271 nm.

JERASH UNIVERSITY
FACULTY OF PHARMACY
MEDICINAL CHEMISTRY I
LAB.VIII

ASSAY OF NALIDIXIC ACID TABLETS

Procedure:-

- 1- Weigh and powder 20 tablets of Nalidixic acid
- 2- To a quantity containing 50 mg of Nalidixic acid add 75 ml of 1M NaOH
- 3- Shake for 3 min.
- 3- Dilute to 100 ml with 1M NaOH, mix, allow to stand for 15 min.
- 4- Take 1 ml and dilute to 100 ml with water
- 5- Measure the absorbance at 334 nm taking 0.01M NaOH as a reference

Calculate the content of $C_{12}H_{12}N_2O_3$ as the value of A (1%, 1cm) at 334 nm.

JERASH UNIVERSITY
FACULTY OF PHARMACY
MEDICINAL CHEMISTRY I
LAB.IX

ASSAY OF DIMENHYDRINATE TABLETS

Procedure:-

- 1- Weigh and powder 20 tablets
- 2- Dissolve a quantity of the powder containing 0.1 gm of Dimenhydrinate as completely as possible in 20 ml of water
- 3- Add 10 ml of 5M Ammonia, mix, extract with successive quantities of 15,15,15,10 and 10 ml of ether and wash with the combined extracts with 10 ml water
- 4- Evaporate the ether, warm the residue with 10 ml of Ethanol 96%, until dissolved, cool, add 50 ml of 0.01M HCl VS
- 5- Titrate the excess of acid with 0.01M NaOH VS using methyl red as indicator.

Each ml of 0.01M HCl is equivalent to 4.7 mg of $C_{17}H_{12}NO, C_7H_7Cl N_4O_2$

JERASH UNIVERSITY
FACULTY OF PHARMACY
MEDICINAL CHEMISTRY I
LAB.X

ASSAY OF CHLORPROPAMIDE TABLETS

Procedure:-

- 1- Weigh and powder 20 tablets
- 2- Shake a quantity of the powder containing 0.25 gm of Chlorpropamide with 40 ml Methanol for 20 min.
- 3- Add sufficient Methanol to produce 50 ml, mix, filter and dilute 5 ml of the filtrate to 100 ml with 0.1M HCl.
- 4- Dilute 4 ml of this solution to 100 ml with 0.1M HCl
- 5- Measure the absorbance of the resulting solution at the maximum at 232 nm

Calculate the content of $C_{10}H_{13}ClN_2O_3S$, taking 598 as the value of $E_1^{1\%}$.

JERASH UNIVERSITY
FACULTY OF PHARMACY
MEDICINAL CHEMISTRY I
LAB.XI

Assay of Propranolol Tablets

- 1- Weight and powder 20 tablets.
- 2- Shake a quantity of powder containing 20 mg of Propranolol Hydrochloride with 20 ml of water for 10 minutes.
- 3- Add 50 ml of methanol, shake for a further 10 minutes, add sufficient methanol to produce 100 ml and filter.
- 4- Dilute 10 ml of the filtrate to 50 ml with methanol
- 5- Measure the absorbance of the resulting solution at the maximum at 290 nm.
- 6- Calculate the content of $C_{16}H_{21}NO_2, HCl$ taking 206 as the value of A (1% , 1 cm) at the maximum at 290 nm.